


Acids

Arrhenius acids


• produce hydrogen ions (H⁺) ions when they dissolve in water

$$HCl(g) \xrightarrow{H_2O(l)} H^+(aq) + Cl^-(aq)$$

- are also electrolytes because they produce H⁺ in water
- have a sour taste
- turn blue litmus red
- corrode some metals














The pH Scale





The pH of a solution

- is used to indicate the acidity of a solution
- has values that usually range from 0 to 14
- is acidic when the values are less than 7
- is neutral at a value of 7
- is **basic** when the values are greater than 7



$$\begin{array}{lll} \text{Acidic solution} & pH < 7.0 & \left[H_3 O^+ \right] > 1 \times 10^{-7} \, \text{M} \\ \text{Neutral solution} & pH = 7.0 & \left[H_3 O^+ \right] = 1 \times 10^{-7} \, \text{M} \\ \text{Basic solution} & pH > 7.0 & \left[H_3 O^+ \right] < 1 \times 10^{-7} \, \text{M} \\ \end{array}$$

$$[H^+] = 10^{-pH}$$

$$pH = -log[H^+]$$







The pH Scale






рН	Hydronium ion concentration	(moles/L)
1	.1	(1 × 10 ⁻¹)
2	.01	(1 × 10 ⁻²)
3	.001	(1 × 10 ⁻³)
4	.0001	(1 × 10 ⁻⁴)
5	.00001	(1 × 10 ⁻⁵)
6	.000001	(1 × 10 ⁻⁶)
7	.0000001	(1 × 10 ⁻⁷)
8	.00000001	(1 × 10 ⁻⁸)
9	.000000001	(1 × 10 ⁻⁹)
10	.000000001 (1 × 10 ⁻¹⁰)
11	.0000000001 (1 × 10 ⁻¹¹)
12	.0000000001 (1 × 10 ⁻¹²)
13	.000000000001 (1 × 10 ⁻¹³)
14	.000000000001 (1 × 10 ⁻¹⁴)





$$[\mathrm{H}^+] = 10^{-\mathrm{pH}}$$

$$pH = -log[H^+]$$





Solutions for A/B

- ▶ 0.1 M NaCl
- ▶ 0.1 M HCl
- 0.1 M acetic acid
- 0.1 M NaOH
- ▶ 0.1 M NH₃

Solutions for C (5 ml each)

- Deionized Water
- 0.1 M NaCl
- ▶ 0.1 M buffer (high pH)
- ▶ 0.1 M buffer (low pH)

Part A: pH Indicators

Use pH paper to estimate the pH of each solution.

Part B: Measuring pH

Use the pH meter to measure and record the pH of each solution.

Part C: pH Buffers

Observe what happens to the pH of each solution when you add small amounts of acid and base.




Solutions for A/B

- ▶ 0.1 M NaCl
- ▶ 0.1 M HCl
- 0.1 M acetic acid
- ▶ 0.1 M NaOH
- ▶ 0.1 M NH₃

Part A: pH Indicators

Use pH paper to estimate the pH of each solution.





Solutions for A/B

- ▶ 0.1 M NaCl
- ▶ 0.1 M HCl
- 0.1 M acetic acid
- ▶ 0.1 M NaOH
- ▶ 0.1 M NH₃

Part B: Measuring pH

Use the pH meter to measure and record the pH of each solution.





Solutions for A/B

- 0.1 M NaCl
- ▶ 0.1 M HCl
- 0.1 M acetic acid
- ▶ 0.1 M NaOH
- ▶ 0.1 M NH₃

Solutions for C (5 ml each)

- 1. Deionized Water
- 2. 0.1 M NaCl
- 3. 0.1 M buffer (high pH)
- 4. 0.1 M buffer (low pH)

Part C: pH Buffers

Observe what happens to the pH of each solution when you add small amounts of acid and base.

ACID TEST

- Prepare 4 test tubes
 - put 5 mL of each solution into a separate tube
- Estimate the pH with pH paper
- Measure the pH with a pH meter
- Add 3 drops of 0.1 M HCl to each tube
- Estimate the pH with pH paper
- Measure the pH with a pH meter

What is the [H+] concentration of each?



Solutions for A/B

- 0.1 M NaCl
- ▶ 0.1 M HCl
- 0.1 M acetic acid
- ▶ 0.1 M NaOH
- ▶ 0.1 M NH₃

Solutions for C (5 ml each)

- 1. Deionized Water
- 2. 0.1 M NaCl
- 3. 0.1 M buffer (high pH)
- 4. 0.1 M buffer (low pH)

Part C: pH Buffers

Observe what happens to the pH of each solution when you add small amounts of acid and base.

BASE TEST

- Prepare 4 test tubes
 - put 5 mL of each solution into a separate tube
- Estimate the pH with pH paper
- Measure the pH with a pH meter
- Add 3 drops of 0.1 M NaOH to each tube
- Estimate the pH with pH paper
- Measure the pH with a pH meter

What is the [H+] concentration of each?



Questions?

